UNISTATION™

2nd Version

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UNISTATION™ 2nd Version

U.S. Department of Health & Human Services

FDA U.S. Food and Drug Administration
Protecting and Promoting Your Health

Home Food Drugs Medical Devices Radiation-Emitting Products Vaccines, Blood & Biologics Animal & Veterinary Cosmetics Tobacco Products

Establishment Registration & Device Listing

FDA Home Medical Devices Databases

New Search Back To Search Results

Proprietary Name:	CENTRicyte 1000 and Accessories
Classification Name:	CENTRIFUGES (MICRO, ULTRA, REFRIGERATED) FOR CLINICAL USE
Product Code:	J0C
Device Class:	1
Regulation Number:	882.2050
Medical Specialty:	Clinical Chemistry
Registered Establishment Name:	NEOGENESIS
Owner/Operator:	NeoGenesis
Owner/Operator Number:	10048015
Establishment Operations:	Specification Developer

Page Last Updated: 12/01/2014
Note: If you need help accessing information in different file formats, see [Instructions for Downloading Viewers and Players](#).

FDA Accessibility Contact FDA Careers FDA Basics FOIA No Fear Act Site Map Transparency Website Policies

U.S. Food and Drug Administration
10903 New Hampshire Avenue
Silver Spring, MD 20993
Ph: 1-888-INFO-FDA (1-888-463-6332)
Email FDA

U.S. Department of Health & Human Services

Combination Products
Advisory Committees
Science & Research
Regulatory Information
Safety
Emergency Preparedness
International Programs
News & Events
Training and Continuing Education
Inspections/Compliance
State & Local Officials
Consumers
Industry
Health Professionals
FDA Archive

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Youtube Link of UNISTATION Introduction



<http://youtu.be/-Jy3t4y1ny0>

Youku Link of UNISTATION Introduction (Chinese)

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

Breast Fat Transplant Breast enlargement with autologous fat

Fat is harvested from suitable localities (bottom, hips) with water-jet liposuction, which is a particularly gentle method with minimal postoperative pain and bruises.

300 mls of pure fat must be collected, at least. Harvested fat is divided into two parts, and then, in local anaesthetics, it is injected into the base of both breasts.

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TISSUE ENGINEERING
Volume 7, Number 2, 2001
Mary Ann Liebert, Inc.

Multilineage Cells from Human Adipose Tissue: Implications for Cell-Based Therapies

PATRICIA A. ZUK, Ph.D.,^{1,2} MIN ZHU, M.D.,^{1,2} HIROSHI MIZUNO, M.D.,²
JERRY HUANG, B.S.,² J. WILLIAM FUTRELL, M.D.,³ ADAM J. KATZ, M.D.,²
PROSPER BENHAIM, M.D.,² H. PETER LORENZ, M.D.,²
and MARC H. HEDRICK, M.D.²

ABSTRACT

Future cell-based therapies such as tissue engineering will benefit from a source of autologous pluripotent stem cells. For mesodermal tissue engineering, one such source of cells is the bone marrow stroma. The bone marrow compartment contains several cell populations, including mesenchymal stem cells (MSCs) that are capable of differentiating into adipogenic, osteogenic, chondrogenic, and myogenic cells. However, autologous bone marrow procurement has potential limitations. An alternate source of autologous adult stem cells that is obtainable in large quantities, under local anesthesia, with minimal discomfort would be advantageous. In this study, we determined if a population of stem cells could be isolated from human adipose tissue. Human adipose tissue, obtained by suction-assisted liposuction (*i.e.*, liposuction), was processed to obtain a fibroblast-like population of cells or a processed lipoaspirate (PLA). These PLA cells can be maintained *in vitro* for extended periods with stable population doubling and low levels of senescence. Immunofluorescence and flow cytometry show that the majority of PLA cells are of mesodermal or mesenchymal origin with low levels of contaminating pericytes, endothelial cells, and smooth muscle cells. Finally, PLA cells differentiate *in vitro* into adipogenic, chondrogenic, myogenic, and osteogenic cells in the presence of lineage-specific induction factors. In conclusion, the data support the hypothesis that a human lipoaspirate contains multipotent cells and may represent an alternative stem cell source to bone marrow-derived MSCs.

INTRODUCTION

THE THERAPEUTIC POTENTIAL of multilineage stem cells for applications such as tissue engineering and gene therapy is enormous. Conceptually, there are two general types of stem cells potentially useful for these applications: embryonic stem cells (ESCs) and autologous stem cells. Although theoretically appealing because of their pluripotentiality, the practical use of ESCs is limited due to potential problems of

Isolation and culture of stem cells—PLA and MSCs

Human adipose tissue was obtained from elective liposuction procedures under local anesthesia (HSPC #98-08 011-02). In this procedure, a hollow blunt-tipped cannula was introduced into the subcutaneous space through small (~1 cm) incisions. The cannula was attached to gentle suction and moved through the adipose compartment, mechanically disrupting the fat tissue. A solution of saline and the vasoconstrictor epinephrine was infused into the adipose compartment to minimize blood loss and contamination of the tissue by peripheral blood cells. The raw lipoaspirate (~300 cc) was processed according to established methodologies to obtain a stromal vascular fraction (SVF).^{33,34} To isolate the SVF, lipoaspirates were washed extensively with equal volumes of phosphate-buffered saline (PBS), and the ECM was digested at 37°C for 30 min with 0.075% collagenase. Enzyme activity was neutralized with Dulbecco's modified Eagle's medium (DMEM), containing 10% FBS and centrifuged at 1200 × *g* for 10 min to obtain a high-density

SVF pellet. The pellet was resuspended in 160 mM NH₄Cl and incubated at room temperature for 10 min to lyse contaminating red blood cells. The SVF was collected by centrifugation, as detailed above, filtered through a 100-μm nylon mesh to remove cellular debris and incubated overnight at 37°C/5% CO₂ in control medium (DMEM, 10% FBS, 1% antibiotic/antimycotic solution). Following incubation, the plates were washed extensively with PBS to remove residual nonadherent red blood cells. The resulting cell population was termed a processed lipoaspirate (PLA), to distinguish it from the SVF obtained from excised adipose tissue. PLA cells were maintained at 37°C/5% CO₂ in noninductive control medium. Cells did not require specific FBS sera lots for expansion and differentiation (data not shown). For immunofluorescence studies, a population of MSCs was obtained from human bone marrow aspirates according to the protocol of Rickard *et al.*¹⁷ and maintained in control medium. To prevent spontaneous differentiation, cells were maintained at subconfluent levels.

¹Dr. Zuk and Dr. Zhu are co-first authors.

²Laboratory for Regenerative Bioengineering and Repair, Departments of Surgery and Orthopaedic Surgery, UCLA School of Medicine, Los Angeles, California.

³Division of Plastic and Reconstructive Surgery, University of Pittsburgh School of Medicine, Pittsburgh, Pennsylvania.

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Reference: Aesth Plast Surg 2008;32:48-55 Dermatol Surg 2008;34:1178-1185

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ORIGINAL ARTICLE

Cell-Assisted Lipotransfer for Cosmetic Breast Augmentation: Supportive Use of Adipose-Derived Stem/Stromal Cells

Kotaro Yoshimura · Katsujiro Sato · Noriyuki Aoi · Masakazu Kurita ·
Toshitsugu Hirohi · Kiyonori Harii

Received: 27 April 2007 / Accepted: 22 May 2007 / Published online: 1 September 2007
© Springer Science+Business Media, LLC 2007

Abstract

Background Lipoinjection is a promising treatment but has some problems, such as unpredictability and a low rate of graft survival due to partial necrosis.

Methods To overcome the problems with lipoinjection, the authors developed a novel strategy known as cell-assisted lipotransfer (CAL). In CAL, autologous adipose-derived stem (stromal) cells (ASCs) are used in combination with lipoinjection. A stromal vascular fraction (SVF) containing ASCs is freshly isolated from half of the aspirated fat and recombined with the other half. This process converts relatively ASC-poor aspirated fat to ASC-rich fat. This report presents the findings for 40 patients who underwent CAL for cosmetic breast augmentation.


Results Final breast volume showed augmentation by 100 to 200 ml after a mean fat amount of 270 ml was injected. Postoperative atrophy of injected fat was minimal and did not change substantially after 2 months. Cyst formation or microcalcification was detected in four patients. Almost all the patients were satisfied with the soft and natural-appearing augmentation.

Conclusions The preliminary results suggest that CAL is effective and safe for soft tissue augmentation and superior to conventional lipoinjection. Additional study is necessary to evaluate the efficacy of this technique further.


Keywords Adipose-derived stem/stromal cells · Breast augmentation · Cell therapy · Lipoinjection · Tissue engineering

*Mean Injection: 270cc
Result after 2 years: Max 200cc
Survival Rate: Max 74%*

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Fig. 4 Clinical views of a patient in group A (patient 1). Preoperative (*left*) and postoperative (*right*) views at 24 months. A 22-year-old woman underwent breast augmentation using cell-assisted lipotransfer (CAL) (290 ml in each breast), with satisfactory results at 24 months. Her breast circumference increased by 5 cm, and her augmented breast mounds remained soft and natural appearing without injection scars or subcutaneous indurations



Fig. 6 Clinical views of a patient in group C (patient 2). Preoperative (*top*) and postoperative (*bottom*) views at 12 months. A 32-year-old woman underwent breast augmentation with cell-assisted lipotransfer (CAL) (280 ml in each breast). Her breast circumference difference increased from 9 cm (baseline) to 14.5 cm (at 12 months). The breast mounds are soft and natural appearing with no visible injection scars


Fig. 8 Clinical views of a patient in group C (patient 3). Preoperative (*top*) and postoperative (*bottom*) views at 24 months. A 30-year-old woman underwent breast augmentation with cell-assisted lipotransfer (CAL) (310 ml in each breast). Her breasts were dramatically augmented with an increase in breast circumference difference by 8 cm at 24 months. The breast mounds were soft with no subcutaneous indurations. An original inframammary fold on the left breast is slightly visible, but injection scars are not visible

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Cell-Assisted Lipotransfer for Facial Lipoatrophy: Efficacy of Clinical Use of Adipose-Derived Stem Cells

KOTARO YOSHIMURA, MD,* KATSUJIRO SATO, MD,[†] NORIYUKI AOI, MD,* MASAKAZU KURITA, MD,[‡] KENTA INOUE, MD,* HIROTAKA SUGA, MD,* HITOMI ETO, MD,* HARUNOSUKE KATO, MD,* TOSHITSUGU HIROH, MD,[§] AND KIYONORI HARIU, MD[‡]

BACKGROUND Lipoinjection is a promising treatment, but its efficacy in recontouring facial lipoatrophy remains to be established.

OBJECTIVE The objective was to evaluate the efficacy and adverse effects of lipoinjection and supplementation of adipose-derived stem/stromal cells (ASCs) to adipose grafts.

METHODS To overcome drawbacks of autologous lipoinjection, we have developed a novel strategy called cell-assisted lipotransfer (CAL). In CAL, stromal vascular fraction containing ASCs was freshly isolated from half of an aspirated fat sample and attached to the other half of aspirated fat sample with the fat acting as a scaffold. This process converts relatively ASC-poor aspirated fat into ASC-rich fat. We performed conventional lipoinjection (non-CAL; $n=3$) or CAL ($n=3$) on six patients with facial lipoatrophy due to lupus profundus or Parry-Romberg syndrome.

RESULTS All patients obtained improvement in facial contour, but the CAL group had a better clinical improvement score than did the non-CAL patients, although the difference did not reach statistical significance ($p=.11$). Adipose necrosis was found in one non-CAL case who took perioperative oral corticosteroids.

CONCLUSION Our results suggest that CAL is both effective and safe and potentially superior to conventional lipoinjection for facial recontouring.

The authors have indicated no significant interest with commercial supporters.

Facial lipoatrophy is a disfiguring and socially disabling problem that accompanies several inherited and acquired diseases (see Ascher et al.¹ for review). Lupus erythematosus profundus and scleroderma en coup de sabre (en coup de sabre morphea) frequently give rise to facial lipoatrophy, which is often the most problematic manifestation of the disease for the patients even though they require medication, such as oral corticosteroids, to suppress other symptoms associated with the systemic morbidities. Facial lipoatrophy is also seen in patients with human immunodeficiency virus (HIV) infection and Parry-Romberg syndrome (hemifacial progressive atrophy or idiopathic hemifacial atrophy), which involves not only facial

lipoatrophy but also bony defects.² Thus far there is no medical treatment to correct facial lipoatrophy.

Lipoinjection has been used for treating lipoatrophy such as that associated with HIV infection with moderate success,^{3,4} while microsurgical tissue transfer has been a standard surgical treatment for recontouring facial defects seen in Parry-Romberg syndrome. However, the surgical procedure leaves conspicuous scars on the face and the donor site.^{5,6} Autologous lipoinjection is a promising treatment for soft tissue augmentation because there is no associated incision scar or complications associated with foreign materials. Although

TABLE 1. Summary of Patients' Data

	Method					
	Non-CAL			CAL		
Case	1	2	3	1	2	3
Sex	Male	Female	Female	Male	Female	Female
Age (year)	25	42	42	35	33	48
Diagnosis	LEP	LEP	LEP	PRS	LEP	LEP
Severity grading	3	5	2	4	5	4
Perioperative oral corticosteroids	—	+	—	—	+	—
Donor site	Abdomen	Abdomen	Abdomen	Abdomen	Abdomen	Abdomen
Volume of injection (mL)	100	250	50	110	90	100
Follow-up periods (month)	10	9	9	13	9	9
Clinical results	Fair	Good	Fair	Excellent	Good	Good
Complications	—	Adipose necrosis	—	—	—	—

LEP, lupus erythematosus profundus; PRS, Parry-Romberg syndrome.

*Department of Plastic Surgery, University of Tokyo School of Medicine, Tokyo; [†]Cellport Clinic Yokohama, Yokohama; [‡]Department of Plastic Surgery, Kyorin University School of Medicine, Tokyo; [§]Ritz Cosmetic Surgery Clinic Tokyo, Tokyo, Japan

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☐ **Semi-automated extraction and characterization of Stromal Vascular Fraction using a new medical device**

Authors: Hanke, Alexander | Pranti, Lukas | Wenzel, Carina | Nerlich, Michael | Brockhoff, Gero | Loibl, Markus | Gehmert, Sebastian

Article Type: Research Article

Abstract: INTRODUCTION: The stem cell rich Stromal Vascular Fraction (SVF) can be harvested by processing lipo-aspirate or fat tissue with an enzymatic digestion followed by centrifugation. To date neither a standardised extraction method for SVF nor a generally admitted protocol for cell application in patients exists. A novel commercially available semi-automated device for the extraction of SVF promises sterility, consistent results and usability in the clinical routine. The aim of this work was to compare the quantity and quality of the SVF between the new system and an established manual laboratory method. MATERIAL AND METHODS: SVF was extracted from lipo-aspirate both ... [Show more](#)

Keywords: Stromal vascular fraction, lipo-aspirate, adipose tissue, mesenchymal stem cells

DOI: 10.3233/CH-168124

Citation: Clinical Hemorheology and Microcirculation, vol. 64, no. 3, pp. 403-412, 2016

- (1) P. Zuk, et al. Multilineage Cells from Human Adipose Tissue: Implications for Cell-Based Therapies. *TISSUE ENGINEERING* 2001;7:212-228.
- (2) I. Hong, et al. Isolation and culture methods of adipose-derived stem cells. *IFATS* 2009 Oral Abstracts;17.
- (3) K. Yoshimura, et al. Cell-Assisted Lipotransfer for Cosmetic Breast Augmentation: Supportive Use of Adipose-Derived Stem/Stromal Cells. *Aesth Plast Surg* 2008;32:48-55.
- (4) K. Yoshimura, et al. Cell-Assisted Lipotransfer for Facial Lipoatrophy: Efficacy of Clinical Use of Adipose-Derived Stem Cells. *Dermatol Surg* 2008;34:1178–1185.
- (5) H. Suga, et al. Numerical Measurement of Viable and Nonviable Adipocytes and Other Cellular Components in Aspirated Fat Tissue. *Plast. Reconstr. Surg.* 2008;122:113.
- (6) J. Mitchell, et al. Immunophenotype of Human Adipose-Derived Cells: Temporal Changes in Stromal-Associated and Stem Cell-Associated Markers. *STEM CELLS* 2006;24:380.
- (7) A Hanke et al. Semi-automated extraction and characterization of Stromal Vascular Fraction using a new medical device. *Clin Hemorheol Microcirc* 64 (3), 403-412. 2016.

NeoGenesis Co., Ltd.
#302, E&C Dream Tower, 146, Seonyu-ro, Yeongdeungpo-gu,
Seoul, South Korea (07255)
Tel. +82-2-583-1348 Fax. +82-2-2038-0848
Homepage: www.neogenesis.co.kr
Email: info@neogenesis.co.kr